

Art.	Descrizione	Description
	<p><b>BILANCE ELETTRONICHE ECONOMICHE SERIE TKW.</b>            Bilance elettroniche di precisione a celle di carico.            Tara automatica su tutta la portata.            Preselezione del limite di peso alto o basso, con allarme acustico.            Memorizzazione pesate effettuate.            Funzione contapezzi.            Unità di misura selezionabile in g-Kg-Oz-lb.            Lettura digitale con ampio display a LCD retroilluminato.            Piatto inox mm. 300x225.            Funzionamento a batterie ricaricabili.            Complete di alimentatore di rete V 230 50/60 Hz.            Dimensioni di ingombro mm. 320x340x125.            Peso netto Kg 3,9.</p> <p><b>Modello TKW-6</b>            Portata 6000 gr., lettura 0,1 g.</p> <p><b>Modello TKW-15++</b>            Portata 15000 gr., lettura 0,1 g.</p> <p><b>Modello TKW-15</b>            Portata 15 Kg, lettura 0,2 g.</p> <p><b>Modello TKW-30</b>            Portata 30 Kg, lettura 0,5 g.</p>	<p><b>LOW COST ELECTRONIC BALANCES, SERIES QHW</b>            Electronic precision balances with load cells.            Automatic tare on the whole capacity.            Pre-selection of max./min. weight thresholds, with acoustic alarm.            Storage of weight readings.            Piece counter function.            Selectable unit of measure in g, Kg, oz or lb.            Wide, backlit digital liquid crystal display.            Stainless steel plate size 300x225 mm.            Powered by rechargeable batteries.            Complete with 230 V, 50/60 Hz feeder.            Overall dimensions: 320x340x125 mm.            Net weight: 3.9 Kg</p> <p><b>Model QHW-6</b>            Capacity 6000 g, reading 0.1 g.</p> <p><b>Model QHW-15++</b>            Capacity 15000 g, reading 0.1 g.</p> <p><b>Model QHW-15</b>            Capacity 15 Kg, reading 0.2 g.</p> <p><b>Model QHW-30</b>            Capacity 30 Kg, reading 0.5 g.</p>



**Serie QHW  
QHW Series**