

Art.	Descrizione	Description
24GB6892/1	<p>CONDUTTIVIMETRO TASCABILE mod.COND 5 Conduttimetro portatile a tenuta stagna IP67 per misura Conducibilità e TDS e Salinità. Cambio scala automatico. Con microprocessore ed elettrodo intercambiabile di elevata durata. Ampio display LCD a 3 colori con indicazione degli standard utilizzati per la taratura e della stabilità della misura. Taratura a pulsante su 3 punti con 3 valori di standard a scelta tra 84, 1413 µS e 12,88 mS. Lettura della temperatura. Compensazione automatica della temperatura con TR= 20/25°C e coefficiente regolabile TC= 0,00...4,00%/°C. Fattore TDS: 0,40...1,00. Messaggi di auto-diagnosi. Alimentazione 4 batterie 1,5V AAA, compatibile anche con batterie ricaricabili. Autospegnimento. Il Kit comprende: 1x55ml soluzioni standard 1413 µS e 12,88 mS, cinghietta, fazzoletti e valigetta.</p> <p>Campo di misura: COND 0...200 µS - 0...2000 µS - 0...20 mS; TDS 0,1ppm...10ppt; Salinità 0,01...10 ppt;</p> <p>Dimensioni: (ø x L) mm.35 x 175. Peso: 130 g</p>	<p>POCKET CONDUCTIVITY-METER mod. COND 5 Conductivity portable meter waterproof IP67 to measure conductivity and TDS and Salinity. Automatic scale change. With microprocessor and interchangeable electrode of high durability. Large LCD display with 3 colours indicating the standard used for calibration and measurement stability. Calibration button on 3 points with 3 values of standard choice between 84, 1413µS and 12.88 mS mS . Temperature reading. Automatic temperature compensation with TR = 20/25 °C and adjustable factor TC = 0.00 ... 4.00% / °C. TDS Factor: 0.40 ... 1.00. Messages of self-diagnosis. Power supply 4 x 1.5V AAA batteries, also compatible with rechargeable batteries. Auto Off. The kit includes: 1x55ml standard solutions in 1413 µS and 12.88 mS, strap, paper tissues, quick manual, 4 x 1.5V AAA batteries and carry case.</p> <p>Measuring range: COND 0 ... 200 µS - 0 ... 2000 µS - 0 ... 20 mS; TDS 0.1ppm ... 10ppt; Salinity 0.01 ... 10 ppt; °C: 0.0 ... 50.0 (0.1 °C)</p> <p>Dimensions: (ø x L) 35 x 175 mm Weight: 130 g</p>
24GB6892/1/A	Sensore di ricambio per COND 5	Spare electrode for COND 5



Art. 24GB6892/1